

# Growth of crystalline C<sub>60</sub> by evaporation

T. Moorsom

This application note describes the growth of crystalline thin films of C<sub>60</sub>.

## 1 Introduction

C<sub>60</sub>, also known as Buckminsterfullerene, is the most common and stable of the fullerenes, comprising a near spherical cage of sixty carbon atoms. C<sub>60</sub> has a wide range of applications in molecular electronics and spintronics, solar cells and molecular magnetism. [1–3] The growth of high quality, crystalline C<sub>60</sub> films is vital for the production of hybrid meta-materials and devices. In this note, we outline the optimal growth parameters for the evaporation of highly crystalline C<sub>60</sub> films from an effusion cell onto Pt substrates.

## 2 Growth

C<sub>60</sub> can be grown easily on various metallic and semi-conducting substrates. C<sub>60</sub> grows poorly on oxides due to high surface tension, causing large clusters to form. [4] Epitaxial platinum with a (111) texture is an ideal metallic substrate for C<sub>60</sub> growth. Pt (111) films were grown on Al<sub>2</sub>O<sub>3</sub> substrates (see application note for e-beam Pt growth). The substrate temperature was maintained at 20° C. C<sub>60</sub> was evaporated from a graphite crucible in a single filament effusion cell. The effusion cell temperature was increased from standby at a rate of 10 °C per minute to 430 °C. The deposition rate was measured to be 0.5 Å/s using a quartz balance. The rate varies during growth, reaching a maximum value of 0.74 Å/s. The substrate was rotated at 90°/s during deposition. the background pressure was  $8 \times 10^{-10}$  mbar. After growth, the film was capped with a 15 nm thick film of Nb to protect it from oxidation.

## 3 Properties

Structural characterisation was obtained using X-ray reflectivity (XRR) and TEM. Figure 1 shows the XRR data taken of the bilayer film. C<sub>60</sub> has a very low density of 1.65 g/cc. However, a structural peak is observable at 11°. Because of the low scattering length of C<sub>60</sub>, this peak is only visible in samples with extremely high crystallinity. A GenX reflectivity fit was performed for this data. The low angle fit does not




Figure 1: X-ray reflectivity data of a platinum C<sub>60</sub> bilayer (red) and GenX fit (black). The structural C<sub>60</sub> peak at 11° is not captured by the low angle fit.

capture the structural peak as this method does not simulate crystal structure.

A lamella,  $80 \pm 10$  Å thick, was cut from the sample using a Dual-Beam Ga Ion FIB. This lamella was then measured using a Titan FEI TEM at 100 kV. The van der Waals lattice of the C<sub>60</sub> is clearly visible with (111) vertical orientation. The interface between the Pt and C<sub>60</sub> is atomically sharp. The lattice spacing of the C<sub>60</sub> film is  $1.05 \pm 0.03$  nm.

The vibrational spectrum of C<sub>60</sub> films was recorded

GenX fitting parameter	Value
Thickness - C <sub>60</sub> (Å)	975±0.7
Density - C <sub>60</sub> (% of bulk)	100±0.01
RMS Roughness - C <sub>60</sub> (Å)	1.8±0.9
Thickness - Pt (Å)	191±0.2
Density - Pt (% of bulk)	100±0.01
RMS Roughness - Pt (Å)	1.7±0.1

Table 1: Structural parameters obtained through the fitting of XRR data with GenX

